

WHY YOUR REGISTRATION LIBRARY  
SHOULD OUTPUT AN ITK TRANSFORM

## LINKS

- Itk wrapper in icon\_registration (lovingly commented)
  - [https://github.com/uncbiag/ICON/blob/master/src/icon\\_registration/itk\\_wrapper.py](https://github.com/uncbiag/ICON/blob/master/src/icon_registration/itk_wrapper.py)
- Notebook to follow along:
  - <https://colab.research.google.com/drive/1NAFFGCD2hh84kfkCQpNhSC4wK4v4MVIq?usp=sharing>

# THE NIGHTMARE SCENARIO

- My registration code takes as input two  $(1, 1, 128, 128, 128)$  tensors  $I_{\text{fixed}}$  and  $I_{\text{moving}}$ , and outputs a  $(1, 3, 128, 128, 128)$  tensor  $\Phi_{\text{inv}}$ .
- $I_{\text{fixed}} \sim \text{torch.nn.functional.grid\_sample}(I_{\text{moving}}, \Phi_{\text{inv}})$ . Success!
- Later, an MD-PHD gives me two DICOM images to register. One is  $130 \times 130 \times 110$  with spacing  $(1, 8, .7333)$  and the other is  $220 \times 220 \times 100$  with spacing  $(.5, .5, .9)$ . Both have orientation  $[[0, 1, 0], [1, 0, 0], [0, 0, -1]]$
- I say, “easy,” and resample the images to  $128 \times 128 \times 128$ , register them, and email the MD-PHD a pickle of a shape  $(1, 3, 128, 128, 128)$  tensor.
- They stare at me like I’m a cat presenting its owner with a dead mouse.

REGISTERED!



## TWO PROBLEMS

- Downstream medical application code needs to correctly handle image spacing, image orientation, and unusual image origins correctly— in particular, landmarks, volume and size measurements, and mesh data are always in physical coordinates
- If we try to manually handle this metadata by, for each downstream task, manually scaling, transposing, and shifting deformation fields, we will make mistakes and spend a lot of time.

This is notably different from segmentations, where it is trivial to just rescale the output labelmap, copy over the metadata, and call it a day

## TWO PRIORITIES

- 1: No matter what, record transforms in a way that allows transforming physical coordinates
- 2: Maybe we can all record transforms in the same way?

- We need a lingua franca datatype for transforms like itk.Image is for images
  - That understands orientation, spacing, origin
  - That has a correct way to compose transforms, warp images, and transport points

# PROPOSAL

- Write a wrapper for each of our registration libraries that can take input images as `itk.Images` with associated metadata, register them, and return an `itk.Transform`.



# ITK.TRANSFORM

- C++ class with various subclasses, wrapped in python
- Composition, warping points, warping images, computing jacobians, inversion of some types of transforms all implemented in physical coordinates, respecting image metadata
- Serializable
- Integrated with ITK ecosystem – Slicer!

## ICON\_REGISTRATION NOW FOLLOWS THIS CONVENTION

- `pip install icon_registration`
- Pretrained model available for knee mri and brain mri
- Following code from icon\_registration
  - [https://github.com/uncbiag/ICON/blob/master/src/icon\\_registration/itk\\_wrapper.py](https://github.com/uncbiag/ICON/blob/master/src/icon_registration/itk_wrapper.py)

```
[ ] import icon_registration.itk_wrapper as itk_wrapper
import icon_registration.pretrained_models as pretrained_models
import itk
import matplotlib.pyplot as plt

[ ] model = pretrained_models.OAI_knees_registration_model()

#Feel free to experiment with different images here
image_A = itk.imread("9487462_20081003_SAG_3D_DESS_RIGHT_11495603_image.nii.gz")
image_B = itk.imread("9225063_20090413_SAG_3D_DESS_RIGHT_12784112_image.nii.gz")

phi_AB, phi_BA = itk_wrapper.register_pair(model, image_A, image_B)
```

Downloading pretrained model (1.2 GB)

```
[ ] interpolator = itk.LinearInterpolateImageFunction.New(image_A)
warped_image_A = itk.resample_image_filter(image_A,
    transform=phi_AB,
    interpolator=interpolator,
    size=itk.size(image_B),
    output_spacing=itk.spacing(image_B),
    output_direction=image_B.GetDirection(),
    output_origin=image_B.GetOrigin()
)
plt.imshow(image_A[80])
plt.show()
plt.imshow(image_B[80])
plt.show()

plt.imshow(warped_image_A[80])
plt.show()

plt.imshow(itk.checker_board_image_filter(warped_image_A, image_B)[80])
plt.show()
```

USER CODE: LOW RISK  
OF BUGS

## CRIMINALLY UNDERRATED ITK FUNCTIONS

These are super useful, whether or not you follow this proposal:

- [Image.TransformPhysicalPointToContinuousIndex](#)
- Image.TransformContinuousIndexToPhysicalPoint

With this proposal, this would become equally useful:

- Transform.TransformPoint

# COMPLEXITY MOVED TO LIBRARY CODE

```
def register_pair(model, image_A, image_B):

    assert( isinstance(image_A, itk.Image))
    assert( isinstance(image_B, itk.Image))
    icon_registration.network_wrappers.adjust_batch_size(model, 1)
    model.to(config.device)

    A_npy = np.array(image_A)
    B_npy = np.array(image_B)
    A_trch = torch.Tensor(A_npy).to(config.device)[None, None]
    B_trch = torch.Tensor(B_npy).to(config.device)[None, None]

    shape = model.identityMap.shape

    A_resized = F.interpolate(A_trch, size=shape[2:], mode="trilinear", align_corners=False)
    B_resized = F.interpolate(B_trch, size=shape[2:], mode="trilinear", align_corners=False)

    with torch.no_grad():
        x = model(A_resized, B_resized)

    phi_AB = model.phi_AB(model.identityMap)[0].cpu()
    phi_BA = model.phi_BA(model.identityMap)[0].cpu()

    return (
        create_itk_transform(phi_AB, model.identityMap, image_A, image_B),
        create_itk_transform(phi_BA, model.identityMap, image_B, image_A)
    )
```

# COMPLEXITY MOVED TO LIBRARY CODE

```
def create_itk_transform(phi, ident, image_A, image_B):

    disp = phi - ident[0].cpu()

    network_shape_list = list(ident.shape[2:])

    dimension = len(network_shape_list)

    scale = torch.Tensor(network_shape_list)
    for _ in network_shape_list:
        scale = scale[:, None]
    disp *= scale
    tr = itk.DisplacementFieldTransform[itk.D, dimension].New()

    disp_itk_format = disp.double().numpy()[list(reversed(range(dimension)))].transpose(list(range(1, dimension + 1)) + [0])

    itk_disp_field = array_to_vector_image(disp_itk_format)

    tr.SetDisplacementField(itk_disp_field)

    to_aligned = resampling_transform(image_A, list(reversed(network_shape_list)))

    from_aligned = resampling_transform(image_B, list(reversed(network_shape_list))).GetInverseTransform()

    phi_AB_itk = itk.CompositeTransform[itk.D, dimension].New()

    phi_AB_itk.PrependTransform(from_aligned)
    phi_AB_itk.PrependTransform(tr)
    phi_AB_itk.PrependTransform(to_aligned)

    return phi_AB_itk
```

## PREFER TO REPRESENT TRANSFORMS AS COMPOSITIONS

- Instead of attempting to post facto modify a deformation field to respect orientation and spacing, it is better to represent a transform as a composition of
  - An affine transform to go from physical coordinates to the scaled coordinates used by our registration algorithm
  - A deformable transform with the deformation field from our neural network
    - This deformable transform must be represented as a vector image with consistent, but not necessarily physical, metadata
  - An affine transform to go from our nonphysical internal coordinates to physical coordinates.
- As a result, it is critical that compositions are a first class citizen of our transform datatype

WORTHWHILE TO WRITE TRICKY LIBRARY  
CODE INSTEAD OF ASKING USERS TO WRITE  
TRICKY APPLICATION CODE!



## SERIALIZATION AND DESERIALIZATION

```
itk.transformwrite([phi_AB], "phi_AB_.hdf5")
```

```
phi_AB_from_disk = itk.transformread("phi_AB_.hdf5")[0]
```

# SHOULD WE USE ITK.TRANSFORM?

- Pros:
  - APIs for warping an image with metadata, composing transforms, and warping physical points are unambiguous
  - Can switch between deformable transform, affine transform, spline transform without changing application code
  - Integrated closely with itk images, which are the most common way of representing medical images with metadata
  - Does not have to integrate into deep learning code
  - Medical people are (begrudgingly) familiar with ITK
  - Correctness is easy to verify
  - Slicer Compatibility
  - Serialization!
- Cons:
  - Python APIs for warping an image, composing transforms, and warping physical points are clunky [itk type system]
  - Has to be a wrapper: cannot integrate into deep learning code / differentiate through

# WHO'S ON BOARD?

- Full compatibility
  - 3D Slicer
  - ITK Registration
  - icon\_registration
  - SimpleITK
- In Flight
  - ITK\_elastix
  - ANTs: internal representation is itk.Transform?
  - Perhaps MONAI?

Speak up and correct me  
on this slide!